

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### Chicken HSP-70/HSPA9(HeatShock Protein 70) ELISA Kit EBT-ELK4696

Article Name	Chicken HSP-70/HSPA9(HeatShock Protein 70) ELISA Kit
Biozol Catalog Number	EBT-ELK4696
Supplier Catalog Number	ELK4696
Alternative Catalog Number	EBT-ELK4696-96, EBT-ELK4696-48, EBT-ELK4696-96X5
Manufacturer	ELK Biotechnology
Category	Kits/Assays
Species Reactivity	Gallus
Concentration	100 ng/mL
Range	1.57-100 ng/mL
Sensitivity	0.53 ng/mL
UniProt	<a href="#">Q5ZM98</a>
Samples	Serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids

Application Notes

Assay Type: Sandwich. Assay length: 3.5h. Research Area: Tumor immunity, Infection immunity, Developmental science, . Test principle: The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided in this kit has been pre-coated with an antibody specific to Chicken HSP-70/HSPA9. Standards or samples are added to the appropriate microtiter plate wells then with a biotin-conjugated antibody specific to Chicken HSP-70/HSPA9. Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain Chicken HSP-70/HSPA9, biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm 10nm. The concentration of Chicken HSP-70/HSPA9 in the samples is then determined by comparing the OD of the samples to the standard curve