

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Pig D2D(D-Dimer) ELISA Kit EBT-ELK5626

Article Name	Pig D2D(D-Dimer) ELISA Kit
Biozol Catalog Number	EBT-ELK5626
Supplier Catalog Number	ELK5626
Alternative Catalog Number	EBT-ELK5626-96, EBT-ELK5626-48, EBT-ELK5626-96X5
Manufacturer	ELK Biotechnology
Category	Kits/Assays
Species Reactivity	Porcine
Concentration	5000 ng/mL
Range	78.13-5000 ng/mL
Sensitivity	25.4 ng/mL
Samples	plasma

Application Notes

Assay Type: Sandwich. Assay length: 3.5h. Research Area: Hematology, Cardiovascular biology,. Test principle: The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided in this kit has been pre-coated with an antibody specific to Pig D2D. Standards or samples are added to the appropriate microtiter plate wells then with a biotin-conjugated antibody specific to Pig D2D. Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain Pig D2D, biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm 10nm. The concentration of Pig D2D in the samples is then determined by comparing the OD of the samples to the standard curve