

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### Cattle TLR4(Toll Like Receptor 4) ELISA Kit EBT-ELK5795

Article Name	Cattle TLR4(Toll Like Receptor 4) ELISA Kit
Biozol Catalog Number	EBT-ELK5795
Supplier Catalog Number	ELK5795
Alternative Catalog Number	EBT-ELK5795-96, EBT-ELK5795-48, EBT-ELK5795-96X5
Manufacturer	ELK Biotechnology
Category	Kits/Assays
Species Reactivity	Bovine
Concentration	40 ng/mL
Range	0.63-40 ng/mL
Sensitivity	0.243 ng/mL
UniProt	<a href="#">Q9GL65</a>
Samples	serum, plasma, tissue homogenates, cell lysates and other biological fluids

Application Notes

Assay Type: Sandwich. Assay length: 3.5h. Research Area: Signal transduction, Enzyme & Kinase, Tumor immunity, Infection immunity, Pulmonology, . Test principle: The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided in this kit has been pre-coated with an antibody specific to Cattle TLR4. Standards or samples are added to the appropriate microtiter plate wells then with a biotin-conjugated antibody specific to Cattle TLR4. Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain Cattle TLR4, biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm 10nm. The concentration of Cattle TLR4 in the samples is then determined by comparing the OD of the samples to the standard curve