

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### **Rabbit NF-kappaB p105(Nuclear factor NF-kappa-B p105 subunit) ELISA Kit EBT-ELK9180**

Article Name	Rabbit NF-kappaB p105(Nuclear factor NF-kappa-B p105 subunit) ELISA Kit
Biozol Catalog Number	EBT-ELK9180
Supplier Catalog Number	ELK9180
Alternative Catalog Number	EBT-ELK9180-96, EBT-ELK9180-48, EBT-ELK9180-96X5
Manufacturer	ELK Biotechnology
Category	Kits/Assays
Species Reactivity	Rabbit
Concentration	20 ng/mL
Range	0.32-20 ng/mL
Sensitivity	0.113 ng/mL
Samples	tissue homogenates, cell lysates, cell culture supernates and other biological fluids

Application Notes	<p>Assay Type: Sandwich. Assay length: 3.5h. Research Area: Signal transduction,. Test principle: The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided in this kit has been pre-coated with an antibody specific to Rabbit NF-kappaBp105. Standards or samples are added to the appropriate microtiter plate wells then with a biotin-conjugated antibody specific to Rabbit NF-kappaBp105. Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain Rabbit NF-kappaBp105, biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm 10nm. The concentration of Rabbit NF-kappaBp105 in the samples is then determined by comparing the OD of the samples to the standard curve</p>
-------------------	---