

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Presporulation Medium (Powder) USB-P8900

| | |
|----------------------------|---|
| Article Name | Presporulation Medium (Powder) |
| Biozol Catalog Number | USB-P8900 |
| Supplier Catalog Number | P8900 |
| Alternative Catalog Number | USB-P8900-500,USB-P8900-2.5,USB-P8900-10 |
| Manufacturer | US Biological |
| Category | Biochemikalien |
| Product Description | Presporulation Medium is used on strains of yeast cells that do not sporulate efficiently, such as those that carry mutations that cause the cells to grow slowly or ferment poorly, are grown for 1 to 2 days before the cells are transferred to Sporula... |