

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### Human Vacuolar protein sorting-associated protein 37A (VPS37A) ELISA Kit ASC-KTE60048

Article Name	Human Vacuolar protein sorting-associated protein 37A (VPS37A) ELISA Kit
Biozol Catalog Number	ASC-KTE60048
Supplier Catalog Number	KTE60048
Alternative Catalog Number	ASC-KTE60048-48, ASC-KTE60048-96
Manufacturer	Abbkine Scientific
Category	Kits/Assays
Application	ELISA
Species Reactivity	Human
Product Description	This Human Vacuolar protein sorting-associated protein 37A (VPS37A) ELISA Kit employs a two-site sandwich ELISA to quantitate VPS37A....
Range	Please inquire
Sensitivity	Please inquire
Tag	VPS37A
NCBI	<a href="#">137492</a>
UniProt	<a href="#">Q8NEZ2</a>
Samples	Cell culture supernatants Serum Plasma Other biological fluids

## Application Notes

This Human Vacuolar protein sorting-associated protein 37A (VPS37A) ELISA Kit employs a two-site sandwich ELISA to quantitate VPS37A in samples. An antibody specific for VPS37A has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and any VPS37A present is bound by the immobilized antibody. After removing any unbound substances, a biotin-conjugated antibody specific for VPS37A is added to the wells. After washing, Streptavidin conjugated Horseradish Peroxidase (HRP) is added to the wells. Following a wash to remove any unbound avidin-enzyme reagent, a substrate solution is added to the wells and color develops in proportion to the amount of VPS37A bound in the initial step. The color development is stopped and the intensity of the color is measured.