

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Human Uveal autoantigen with coiled-coil domains and ankyrin repeats (UACA) ELISA Kit ASC-KTE60107

Article Name	Human Uveal autoantigen with coiled-coil domains and ankyrin repeats (UACA) ELISA Kit
Biozol Catalog Number	ASC-KTE60107
Supplier Catalog Number	KTE60107
Alternative Catalog Number	ASC-KTE60107-48, ASC-KTE60107-96
Manufacturer	Abbkine Scientific
Category	Kits/Assays
Application	ELISA
Species Reactivity	Human
Product Description	This Human Uveal autoantigen with coiled-coil domains and ankyrin repeats (UACA) ELISA Kit employs a two-site sandwich ELISA to quantitate UACA....
Range	Please inquire
Sensitivity	Please inquire
Tag	UACA
NCBI	55075
UniProt	Q9BZF9
Samples	Cell culture supernatantsSerumPlasmaOther biological fluids

Application Notes	<p>This Human Uveal autoantigen with coiled-coil domains and ankyrin repeats (UACA) ELISA Kit employs a two-site sandwich ELISA to quantitate UACA in samples. An antibody specific for UACA has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and any UACA present is bound by the immobilized antibody. After removing any unbound substances, a biotin-conjugated antibody specific for UACA is added to the wells. After washing, Streptavidin conjugated Horseradish Peroxidase (HRP) is added to the wells. Following a wash to remove any unbound avidin-enzyme reagent, a substrate solution is added to the wells and color develops in proportion to the amount of UACA bound in the initial step. The color development is stopped and the intensity of the color is measured.</p>
-------------------	--