

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### Human Tumor necrosis factor alpha-induced protein 2 (TNFAIP2) ELISA Kit ASC-KTE60235

Article Name	Human Tumor necrosis factor alpha-induced protein 2 (TNFAIP2) ELISA Kit
Biozol Catalog Number	ASC-KTE60235
Supplier Catalog Number	KTE60235
Alternative Catalog Number	ASC-KTE60235-48, ASC-KTE60235-96
Manufacturer	Abbkine Scientific
Category	Kits/Assays
Application	ELISA
Species Reactivity	Human
Product Description	This Human Tumor necrosis factor alpha-induced protein 2 (TNFAIP2) ELISA Kit employs a two-site sandwich ELISA to quantitate TNFAIP2....
Range	Please inquire
Sensitivity	Please inquire
Tag	TNFAIP2
NCBI	<a href="#">7127</a>
UniProt	<a href="#">Q03169</a>
Samples	Cell culture supernatantsSerumPlasmaOther biological fluids

Application Notes

This Human Tumor necrosis factor alpha-induced protein 2 (TNFAIP2) ELISA Kit employs a two-site sandwich ELISA to quantitate TNFAIP2 in samples. An antibody specific for TNFAIP2 has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and any TNFAIP2 present is bound by the immobilized antibody. After removing any unbound substances, a biotin-conjugated antibody specific for TNFAIP2 is added to the wells. After washing, Streptavidin conjugated Horseradish Peroxidase (HRP) is added to the wells. Following a wash to remove any unbound avidin-enzyme reagent, a substrate solution is added to the wells and color develops in proportion to the amount of TNFAIP2 bound in the initial step. The color development is stopped and the intensity of the color is measured.