

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Human Adiponectin (ADP) ELISA Kit ASC-KTE60939

| | |
|----------------------------|--|
| Article Name | Human Adiponectin (ADP) ELISA Kit |
| Biozol Catalog Number | ASC-KTE60939 |
| Supplier Catalog Number | KTE60939 |
| Alternative Catalog Number | ASC-KTE60939-48, ASC-KTE60939-96 |
| Manufacturer | Abbkine Scientific |
| Category | Kits/Assays |
| Application | ELISA |
| Species Reactivity | Human |
| Product Description | Human Adiponectin (ADP) ELISA Kit employs a two-site sandwich ELISA to quantitate Adiponectin in samples.... |
| Range | 125 pg/mL-2000 pg/mL |
| Sensitivity | 10 pg/mL |
| Tag | Adiponectin |
| Samples | Cell culture supernatantsSerumPlasmaOther biological fluids |

Application Notes

Human Adiponectin (ADP) ELISA Kit employs a two-site sandwich ELISA to quantitate Adiponectin in samples. An antibody specific for Adiponectin has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and any Adiponectin present is bound by the immobilized antibody. After removing any unbound substances, HRP-Conjugate Human Adiponectin detection antibody is added to the wells. Following a wash to remove any unbound HRP reagent, a Chromogen solution is added to the wells and color develops in proportion to the amount of Adiponectin bound in the initial step. The color development is stopped and the intensity of the color is measured.