

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Human N-acetyltransferase 6 (NAT6) ELISA Kit ASC-KTE61365

Article Name	Human N-acetyltransferase 6 (NAT6) ELISA Kit
Biozol Catalog Number	ASC-KTE61365
Supplier Catalog Number	KTE61365
Alternative Catalog Number	ASC-KTE61365-48, ASC-KTE61365-96
Manufacturer	Abbkine Scientific
Category	Kits/Assays
Application	ELISA
Species Reactivity	Human
Product Description	This Human N-acetyltransferase 6 (NAT6) ELISA Kit employs a two-site sandwich ELISA to quantitate NAT6....
Range	Please inquire
Sensitivity	Please inquire
Tag	NAT6
NCBI	24142
UniProt	Q93015
Samples	Cell culture supernatantsSerumPlasmaOther biological fluids

Application Notes

This Human N-acetyltransferase 6 (NAT6) ELISA Kit employs a two-site sandwich ELISA to quantitate NAT6 in samples. An antibody specific for NAT6 has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and any NAT6 present is bound by the immobilized antibody. After removing any unbound substances, a biotin-conjugated antibody specific for NAT6 is added to the wells. After washing, Streptavidin conjugated Horseradish Peroxidase (HRP) is added to the wells. Following a wash to remove any unbound avidin-enzyme reagent, a substrate solution is added to the wells and color develops in proportion to the amount of NAT6 bound in the initial step. The color development is stopped and the intensity of the color is measured.