

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### Human NADH-ubiquinone oxidoreductase chain 2 (MT-ND2) ELISA Kit ASC-KTE61490

Article Name	Human NADH-ubiquinone oxidoreductase chain 2 (MT-ND2) ELISA Kit
Biozol Catalog Number	ASC-KTE61490
Supplier Catalog Number	KTE61490
Alternative Catalog Number	ASC-KTE61490-48, ASC-KTE61490-96
Manufacturer	Abbkine Scientific
Category	Kits/Assays
Application	ELISA
Species Reactivity	Human
Product Description	This Human NADH-ubiquinone oxidoreductase chain 2 (MT-ND2) ELISA Kit employs a two-site sandwich ELISA to quantitate MT-ND2....
Range	Please inquire
Sensitivity	Please inquire
Tag	MT-ND2
NCBI	<a href="#">4536</a>
UniProt	<a href="#">P03891</a>
Samples	Cell culture supernatantsSerumPlasmaOther biological fluids

Application Notes

This Human NADH-ubiquinone oxidoreductase chain 2 (MT-ND2) ELISA Kit employs a two-site sandwich ELISA to quantitate MT-ND2 in samples. An antibody specific for MT-ND2 has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and any MT-ND2 present is bound by the immobilized antibody. After removing any unbound substances, a biotin-conjugated antibody specific for MT-ND2 is added to the wells. After washing, Streptavidin conjugated Horseradish Peroxidase (HRP) is added to the wells. Following a wash to remove any unbound avidin-enzyme reagent, a substrate solution is added to the wells and color develops in proportion to the amount of MT-ND2 bound in the initial step. The color development is stopped and the intensity of the color is measured.