

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Mouse Troponin I, fast skeletal muscle (TNNI2) ELISA Kit ASC-KTE70099

Article Name	Mouse Troponin I, fast skeletal muscle (TNNI2) ELISA Kit
Biozol Catalog Number	ASC-KTE70099
Supplier Catalog Number	KTE70099
Alternative Catalog Number	ASC-KTE70099-48, ASC-KTE70099-96
Manufacturer	Abbkine Scientific
Category	Kits/Assays
Application	ELISA
Species Reactivity	Mouse
Product Description	This Mouse Troponin I, fast skeletal muscle (TNNI2) ELISA Kit employs a two-site sandwich ELISA to quantitate TNNI2....
Range	Please inquire
Sensitivity	Please inquire
Tag	TNNI2
NCBI	21953
UniProt	P13412
Samples	Cell culture supernatants Serum Plasma Other biological fluids

Application Notes

This Mouse Troponin I, fast skeletal muscle (TNNI2) ELISA Kit employs a two-site sandwich ELISA to quantitate TNNI2 in samples. An antibody specific for TNNI2 has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and anyTNNI2 present is bound by the immobilized antibody. After removing any unbound substances, a biotin-conjugated antibody specific for TNNI2 is added to the wells. After washing, Streptavidin conjugated Horseradish Peroxidase (HRP) is added to the wells. Following a wash to remove any unbound avidin-enzyme reagent, a substrate solution is added to the wells and color develops in proportion to the amount of TNNI2 bound in the initial step. The color development is stopped and the intensity of the color is measured.