

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Mouse Tropomodulin-2 (TMOD2) ELISA Kit ASC-KTE70155

| | |
|----------------------------|---|
| Article Name | Mouse Tropomodulin-2 (TMOD2) ELISA Kit |
| Biozol Catalog Number | ASC-KTE70155 |
| Supplier Catalog Number | KTE70155 |
| Alternative Catalog Number | ASC-KTE70155-48, ASC-KTE70155-96 |
| Manufacturer | Abbkine Scientific |
| Category | Kits/Assays |
| Application | ELISA |
| Species Reactivity | Mouse |
| Product Description | This Mouse Tropomodulin-2 (TMOD2) ELISA Kit employs a two-site sandwich ELISA to quantitate TMOD2.... |
| Range | Please inquire |
| Sensitivity | Please inquire |
| Tag | TMOD2 |
| NCBI | 50876 |
| UniProt | Q9JKK7 |
| Samples | Cell culture supernatantsSerumPlasmaOther biological fluids |

| | |
|-------------------|--|
| Application Notes | <p>This Mouse Tropomodulin-2 (TMOD2) ELISA Kit employs a two-site sandwich ELISA to quantitate TMOD2 in samples. An antibody specific for TMOD2 has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and any TMOD2 present is bound by the immobilized antibody. After removing any unbound substances, a biotin-conjugated antibody specific for TMOD2 is added to the wells. After washing, Streptavidin conjugated Horseradish Peroxidase (HRP) is added to the wells. Following a wash to remove any unbound avidin-enzyme reagent, a substrate solution is added to the wells and color develops in proportion to the amount of TMOD2 bound in the initial step. The color development is stopped and the intensity of the color is measured.</p> |
|-------------------|--|