

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Mouse Toll-like receptor 5 (TLR-5) ELISA Kit ASC-KTE70209

| | |
|----------------------------|--|
| Article Name | Mouse Toll-like receptor 5 (TLR-5) ELISA Kit |
| Biozol Catalog Number | ASC-KTE70209 |
| Supplier Catalog Number | KTE70209 |
| Alternative Catalog Number | ASC-KTE70209-48, ASC-KTE70209-96 |
| Manufacturer | Abbkine Scientific |
| Category | Kits/Assays |
| Application | ELISA |
| Species Reactivity | Mouse |
| Product Description | This Mouse Toll-like receptor 5 (TLR-5) ELISA Kit employs a two-site sandwich ELISA to quantitate TLR5.... |
| Range | Please inquire |
| Sensitivity | Please inquire |
| Tag | TLR5 |
| NCBI | 53791 |
| UniProt | Q9JLF7 |
| Samples | Cell culture supernatants Serum Plasma Other biological fluids |

Application Notes

This Mouse Toll-like receptor 5 (TLR-5) ELISA Kit employs a two-site sandwich ELISA to quantitate TLR5 in samples. An antibody specific for TLR5 has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and any TLR5 present is bound by the immobilized antibody. After removing any unbound substances, a biotin-conjugated antibody specific for TLR5 is added to the wells. After washing, Streptavidin conjugated Horseradish Peroxidase (HRP) is added to the wells. Following a wash to remove any unbound avidin-enzyme reagent, a substrate solution is added to the wells and color develops in proportion to the amount of TLR5 bound in the initial step. The color development is stopped and the intensity of the color is measured.