

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### Mouse Dehydroepiandrosterone sulfate (DHEA-S) ELISA Kit ASC-KTE71292

Article Name	Mouse Dehydroepiandrosterone sulfate (DHEA-S) ELISA Kit
Biozol Catalog Number	ASC-KTE71292
Supplier Catalog Number	KTE71292
Alternative Catalog Number	ASC-KTE71292-48, ASC-KTE71292-96
Manufacturer	Abbkine Scientific
Category	Kits/Assays
Application	ELISA
Species Reactivity	Mouse
Product Description	This Mouse Dehydroepiandrosterone sulfate (DHEA-S) ELISA Kit employs a two-site sandwich ELISA to quantitate DHEA-S....
Range	Please inquire
Sensitivity	Please inquire
Tag	DHEA-S
Samples	Cell culture supernatants Serum Plasma Other biological fluids

Application Notes

This Mouse Dehydroepiandrosterone sulfate (DHEA-S) ELISA Kit employs a two-site sandwich ELISA to quantitate DHEA-S in samples. An antibody specific for DHEA-S has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and anyDHEA-S present is bound by the immobilized antibody. After removing any unbound substances, a biotin-conjugated antibody specific for DHEA-S is added to the wells. After washing, Streptavidin conjugated Horseradish Peroxidase (HRP) is added to the wells. Following a wash to remove any unbound avidin-enzyme reagent, a substrate solution is added to the wells and color develops in proportion to the amount of DHEA-S bound in the initial step. The color development is stopped and the intensity of the color is measured.