

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### ExKine(TM) Membrane and Cytoplasmic Protein Extraction kit ASC-KTP3005

Article Name	ExKine(TM) Membrane and Cytoplasmic Protein Extraction kit
Biozol Catalog Number	ASC-KTP3005
Supplier Catalog Number	KTP3005
Alternative Catalog Number	ASC-KTP3005-50, ASC-KTP3005-200
Manufacturer	Abbkine Scientific
Category	Kits/Assays
Product Description	Abbkine ExKine(TM) Membrane and Cytoplasmic Protein Extraction Kit enable stepwise separation and preparation of cytoplasmic and membrane protein....
Tag	Cytoplasmic
Application Notes	Abbkine ExKine(TM) Membrane and Cytoplasmic Protein Extraction Kit enable stepwise separation and preparation of non-denatured, active cytoplasmic and membrane protein from mammalian cultured cells or tissues, which can be used directly in a variety of proteomics application. The cells are first broken by homogenization to allow the release of cytosolic and membrane proteins, after which membrane proteins are obtained by high speed centrifugation. Extraction efficiencies and yields will vary depending on cell type as well as the number of times the integral membrane protein spans the lipid bilayer.