

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Histone H3, phosphorylated (Ser10) Sandwich ELISA Kit, BioAssay(TM) USB-H5110-12N

| | |
|----------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Article Name | Histone H3, phosphorylated (Ser10) Sandwich ELISA Kit, BioAssay(TM) |
| Biozol Catalog Number | USB-H5110-12N |
| Supplier Catalog Number | H5110-12N |
| Alternative Catalog Number | USB-H5110-12N-1 |
| Manufacturer | US Biological |
| Category | Kits/Assays |
| Product Description | Detects endogenous levels of Phospho-Histone H3 (Ser10). Using this Sandwich ELISA Kit Phospho- Histone H3 (Ser10) is detected when treated with Calyculin A in NIH 3T3 cells. However, the levels of Histone H3 remains unchanged, as shown by Western an... |