

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Loading Dye Buffer, 6X USB-L3350

Article Name	Loading Dye Buffer, 6X
Biozol Catalog Number	USB-L3350
Supplier Catalog Number	L3350
Alternative Catalog Number	USB-L3350-5X1
Manufacturer	US Biological
Category	Molekularbiologie
Product Description	Loading dye for DNA markers and samples on agarose or polyacrylamide gels. It contains two dyes for easy visual tracking. Recommended Dilutions: Add 1:5 of L3350 to DNA sample, mix 1 volume of dye with 5 volumes of DNA sample. In 1% agarose gels brom...
Form	Supplied as a liquid in 10mM Tris-HCl, pH 7.6, 0.03% bromophenol blue, 0.03% xylene cyanol FF, 60mM EDTA, 60% glycerol.