

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### **MD-2, Control Peptide (MD-2 Protein, ESOP1, ESOP-1, Lymphocyte Antigen 96, LY96) USB-M2767-20D**

|                            |                                                                                                                                                                                                                                                               |
|----------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Article Name               | MD-2, Control Peptide (MD-2 Protein, ESOP1, ESOP-1, Lymphocyte Antigen 96, LY96)                                                                                                                                                                              |
| Biozol Catalog Number      | USB-M2767-20D                                                                                                                                                                                                                                                 |
| Supplier Catalog Number    | M2767-20D                                                                                                                                                                                                                                                     |
| Alternative Catalog Number | USB-M2767-20D-50                                                                                                                                                                                                                                              |
| Manufacturer               | US Biological                                                                                                                                                                                                                                                 |
| Category                   | Molekularbiologie                                                                                                                                                                                                                                             |
| Application                | ELISA, WB                                                                                                                                                                                                                                                     |
| Product Description        | Sequence: 13aa near the middle region of human MD-2 (1). MD-2 enables TLR4 to respond to a wide variety of LPS partial structures, Gram-negative bacteria, and Gram-positive lipoteichoic acid, but not to Gram-positive bacteria, peptidoglycan, and lipo... |
| Purity                     | Purified                                                                                                                                                                                                                                                      |
| Form                       | Supplied as a liquid in PBS, pH 7.2, 0.1% BSA, 0.02% sodium azide.                                                                                                                                                                                            |