

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### **Tissue, Array, Human Tumor, Array II (192 Specimens/47 Pairs) (Paraffin) HisTek(TM) USB-T5595-6335**

|                            |                                                                                                                                                                                                                                                               |
|----------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Article Name               | Tissue, Array, Human Tumor, Array II (192 Specimens/47 Pairs) (Paraffin) HisTek(TM)                                                                                                                                                                           |
| Biozol Catalog Number      | USB-T5595-6335                                                                                                                                                                                                                                                |
| Supplier Catalog Number    | T5595-6335                                                                                                                                                                                                                                                    |
| Alternative Catalog Number | USB-T5595-6335-2,USB-T5595-6335-5                                                                                                                                                                                                                             |
| Manufacturer               | US Biological                                                                                                                                                                                                                                                 |
| Category                   | Sonstiges                                                                                                                                                                                                                                                     |
| Product Description        | Tissues were fixed in formalin after excision, and embedded in paraffin. Paraffin tissue arrays are 5-10um in thickness, and mounted on positively charged glass slides. Each tissue sample is about 2mm in diameter. A reference chart shows all the tiss... |