

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Bovine Annexin V-Fluorescein Apoptosis Assay Kit BYT-ORB1409913

Article Name	Bovine Annexin V-Fluorescein Apoptosis Assay Kit
Biozol Catalog Number	BYT-ORB1409913
Supplier Catalog Number	orb1409913
Alternative Catalog Number	BYT-ORB1409913-500
Manufacturer	Biorbyt
Category	Kits/Assays
Product Description	For quickly and easily distinguishing two populations of dying cells from viable cells using recombinant fluorescein-conjugated Bovine Annexin V and Propidium Iodide. Results can be analyzed via flow cytometer....
Samples	Cell culture
Application Notes	Application Notes: 1. Prepare samples and controls., 2. Dilute 10X Binding Buffer 1:10 WITH DiH2O., 3. Wash cells twice in ice-cold culture medium or PBS and resuspend in ice-cold 1X Binding Buffer at 5×10^5 - 1×10^6 cells/mL., 4. Reconstitute Canine Annexin V-Fluorescein in 1 mL Annexin Reconstitution Buffer, 1X., 5. Dilute reconstituted Annexin V-Fluorescein 1:5 in PBS., 6. Add diluted Annexin V-Fluorescein to each sample at ~1:10., 7. Dilute PI 1:10 in PBS., 8. Add diluted PI to each sample at ~1:20., 9. Incubate 15 minutes on ice and protected from light., 10. Add 400 μ L 1X Binding Buffer to each sample., 11. Analyze with a flow cytometer. Annexin V-Fluorescein excites at 494 nm and emits at 521 nm. PI excites at 536 nm and emits at 617 nm