

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

FAM-Phe-DAP Green FLISP Assay Kit BYT-ORB1409921

Article Name	FAM-Phe-DAP Green FLISP Assay Kit
Biozol Catalog Number	BYT-ORB1409921
Supplier Catalog Number	orb1409921
Alternative Catalog Number	BYT-ORB1409921-100
Manufacturer	Biorbyt
Category	Kits/Assays
Product Description	Quantify intracellular active chymotrypsin-like serine protease activity in vitro. The DAP group phosphorylates the active site serine residue. Analyze your samples using a flow cytometer, fluorescence microscope, or fluorescence plate reader....
Samples	Cell culture, tissue
Application Notes	Application Notes: 1. Add 1 part Antigen-Down HRP Conjugate Stabilizer to 4 parts deionized water and mix for 5-15 minutes., 2. Add the anti-IgG-HRP conjugate directly into the conjugate diluent to a working 1X dilution and mix for 5-15 minutes., 3. Pipette diluted conjugate onto the plate, or store protected from light at 2-8C for use at a later date. 1X working conjugate dilutions can be stored at 2-8C and protected from light, providing a useful shelf-life of 6-12 months