

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

FAM-Leu-DAP Green FLISP Assay Kit BYT-ORB1409928

Article Name	FAM-Leu-DAP Green FLISP Assay Kit
Biozol Catalog Number	BYT-ORB1409928
Supplier Catalog Number	orb1409928
Alternative Catalog Number	BYT-ORB1409928-25
Manufacturer	Biorbyt
Category	Kits/Assays
Product Description	Quantify intracellular active chymotrypsin-like serine protease activity in vitro. The DAP group phosphorylates the active site serine residue. Analyze your samples using a flow cytometer, fluorescence microscope, or fluorescence plate reader....
Samples	Cell culture, tissue
Application Notes	Application Notes: 1. Prepare samples and controls, 2. Dilute cellular wash buffer 1:10 with diH2O, 3. Reconstitute FLISP with 50 µL DMSO, 4. Dilute FLISP 1:5 by adding 200 µL PBS, 5. Add 10 µL FLISP to each sample (~500 µL aliquot of cultured cells), 6. Incubate ~ 1 hour., 7. Remove media and wash cells: add wash buffer and spin cells (twice), or add fresh media and incubate 1 hour, 8. If desired, label with additional stains, such as Hoechst, PI, or an antibody, 9. If desired, fix or embed cells, 10. Analyze with a fluorescence microscope, fluorescence plate reader, or flow cytometer. FAM-FLISP excites at 488 nm and emits at 520 nm