

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Folate Binding Protein Antibody, Unconjugated, Goat, Polyclonal BYT-ORB344182

Article Name	Folate Binding Protein Antibody, Unconjugated, Goat, Polyclonal
Biozol Catalog Number	BYT-ORB344182
Supplier Catalog Number	orb344182
Alternative Catalog Number	BYT-ORB344182-100
Manufacturer	Biorbyt
Host	Goat
Category	Antikörper
Application	ELISA, IP, WB
Species Reactivity	Bovine
Immunogen	Folate Binding Protein [Bovine Milk]
Conjugation	Unconjugated
Product Description	Folate Binding Protein antibody...
Clonality	Polyclonal
Concentration	1 mg/ml
NCBI	001193456
UniProt	P02702
Buffer	Preservative: 0.01% (w/v) Sodium Azide. Stabilizer: None, Buffer: 0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2

Purity	Anti-FOLATE BINDING PROTEIN is an IgG fraction antibody purified from monospecific antiserum by a multi-step process which includes delipidation, salt fractionation and ion exchange chromatography followed by extensive dialysis against the buffer stated above. Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-Goat Serum as well as purified and partially purified Folate Binding Protein [Bovine Milk]. Cross reactivity against Folate Binding Protein from other tissues and species may occur but have not been specifically determined.
Form	Lyophilized
Application Dilute	ELISA: 1:5,000 - 1:20,000, IP: 1:100, WB: 1:500 - 1:5,000
Application Notes	Application Notes: FOLATE BINDING PROTEIN Antibody has been assayed against 1.0 ug of Folate Binding Protein [Bovine Milk] in a standard ELISA using Peroxidase conjugated Affinity Purified anti-Goat IgG [H&L] (Rabbit) and (ABTS (2,2-azino-bis-[3-ethylbenthiazoline-6-sulfonic acid]) as a substrate for 30 minutes at room temperature. A working dilution of 1:10,000 to 1:40,000 of the reconstitution concentration is suggested for this product. Highest levels of antigen binding and antibody reactivity were observed using NuncO brand MaxiSorp immunoplates