

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Rat Mitochondrial import receptor subunit TOM40B (TOMM40L) ELISA Kit BYT-ORB3145658

| | |
|----------------------------|---|
| Article Name | Rat Mitochondrial import receptor subunit TOM40B (TOMM40L) ELISA Kit |
| Biozol Catalog Number | BYT-ORB3145658 |
| Supplier Catalog Number | orb3145658 |
| Alternative Catalog Number | BYT-ORB3145658-48, BYT-ORB3145658-5, BYT-ORB3145658-50, BYT-ORB3145658-96 |
| Manufacturer | Biorbyt |
| Category | Kits/Assays |
| Application | ELISA |
| Species Reactivity | Rat |
| Product Description | This Rat Mitochondrial import receptor subunit TOM40B (TOMM40L) ELISA Kit employs a two-site sandwich ELISA to quantitate TOMM40L.... |
| Range | Please inquire |
| Sensitivity | Please inquire |
| Samples | Cell culture supernatantsSerumPlasmaOther biological fluids |

Application Notes

Application Notes: This Rat Mitochondrial import receptor subunit TOM40B (TOMM40L) ELISA Kit employs a two-site sandwich ELISA to quantitate TOMM40L in samples. An antibody specific for TOMM40L has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and any TOMM40L present is bound by the immobilized antibody. After removing any unbound substances, a biotin-conjugated antibody specific for TOMM40L is added to the wells. After washing, Streptavidin conjugated Horseradish Peroxidase (HRP) is added to the wells. Following a wash to remove any unbound avidin-enzyme reagent, a substrate solution is added to the wells and color develops in proportion to the amount of TOMM40L bound in the initial step. The color development is stopped and the intensity of the color is measured