

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Human Vasodilator-stimulated phosphoprotein (VASP) ELISA Kit BYT-ORB3145798

Article Name	Human Vasodilator-stimulated phosphoprotein (VASP) ELISA Kit
Biozol Catalog Number	BYT-ORB3145798
Supplier Catalog Number	orb3145798
Alternative Catalog Number	BYT-ORB3145798-48,BYT-ORB3145798-96,BYT-ORB3145798-5,BYT-ORB3145798-50
Manufacturer	Biorbyt
Category	Kits/Assays
Application	ELISA
Species Reactivity	Human
Product Description	This Human Vasodilator-stimulated phosphoprotein (VASP) ELISA Kit employs a two-site sandwich ELISA to quantitate VASP....
Range	Please inquire
Sensitivity	Please inquire
Samples	Cell culture supernatantsSerumPlasmaOther biological fluids

Application Notes

Application Notes: This Human Vasodilator-stimulated phosphoprotein (VASP) ELISA Kit employs a two-site sandwich ELISA to quantitate VASP in samples. An antibody specific for VASP has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and any VASP present is bound by the immobilized antibody. After removing any unbound substances, a biotin-conjugated antibody specific for VASP is added to the wells. After washing, Streptavidin conjugated Horseradish Peroxidase (HRP) is added to the wells. Following a wash to remove any unbound avidin-enzyme reagent, a substrate solution is added to the wells and color develops in proportion to the amount of VASP bound in the initial step. The color development is stopped and the intensity of the color is measured