

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### Human Ubiquitin-conjugating enzyme E2 L3 (UBE2L3) ELISA Kit BYT-ORB3145826

|                            |   |
|----------------------------|---|
| Article Name               | Human Ubiquitin-conjugating enzyme E2 L3 (UBE2L3) ELISA Kit   |
| Biozol Catalog Number      | BYT-ORB3145826  |
| Supplier Catalog Number    | orb3145826  |
| Alternative Catalog Number | BYT-ORB3145826-48,BYT-ORB3145826-96,BYT-ORB3145826-5,BYT-ORB3145826-50  |
| Manufacturer               | Biorbyt   |
| Category                   | Kits/Assays   |
| Application                | ELISA   |
| Species Reactivity         | Human   |
| Product Description        | This Human Ubiquitin-conjugating enzyme E2 L3 (UBE2L3) ELISA Kit employs a two-site sandwich ELISA to quantitate UBE2L3.... |
| Range                      | Please inquire  |
| Sensitivity                | Please inquire  |
| Samples                    | Cell culture supernatantsSerumPlasmaOther biological fluids   |

Application Notes

Application Notes: This Human Ubiquitin-conjugating enzyme E2 L3 (UBE2L3) ELISA Kit employs a two-site sandwich ELISA to quantitate UBE2L3 in samples. An antibody specific for UBE2L3 has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and any UBE2L3 present is bound by the immobilized antibody. After removing any unbound substances, a biotin-conjugated antibody specific for UBE2L3 is added to the wells. After washing, Streptavidin conjugated Horseradish Peroxidase (HRP) is added to the wells. Following a wash to remove any unbound avidin-enzyme reagent, a substrate solution is added to the wells and color develops in proportion to the amount of UBE2L3 bound in the initial step. The color development is stopped and the intensity of the color is measured