

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### Mouse TRAF2 and NCK-interacting protein kinase (TNIK) ELISA Kit BYT-ORB3147762

Article Name	Mouse TRAF2 and NCK-interacting protein kinase (TNIK) ELISA Kit
Biozol Catalog Number	BYT-ORB3147762
Supplier Catalog Number	orb3147762
Alternative Catalog Number	BYT-ORB3147762-48, BYT-ORB3147762-5, BYT-ORB3147762-50, BYT-ORB3147762-96
Manufacturer	Biorbyt
Category	Kits/Assays
Application	ELISA
Species Reactivity	Mouse
Product Description	This Mouse TRAF2 and NCK-interacting protein kinase (TNIK) ELISA Kit employs a two-site sandwich ELISA to quantitate TNIK....
Range	Please inquire
Sensitivity	Please inquire
Samples	Cell culture supernatantsSerumPlasmaOther biological fluids

Application Notes

Application Notes: This Mouse TRAF2 and NCK-interacting protein kinase (TNIK) ELISA Kit employs a two-site sandwich ELISA to quantitate TNIK in samples. An antibody specific for TNIK has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and any TNIK present is bound by the immobilized antibody. After removing any unbound substances, a biotin-conjugated antibody specific for TNIK is added to the wells. After washing, Streptavidin conjugated Horseradish Peroxidase (HRP) is added to the wells. Following a wash to remove any unbound avidin-enzyme reagent, a substrate solution is added to the wells and color develops in proportion to the amount of TNIK bound in the initial step. The color development is stopped and the intensity of the color is measured