

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

General Tyrosol (Tyr) ELISA Kit BYT-ORB782099

| | |
|----------------------------|---|
| Article Name | General Tyrosol (Tyr) ELISA Kit |
| Biozol Catalog Number | BYT-ORB782099 |
| Supplier Catalog Number | orb782099 |
| Alternative Catalog Number | BYT-ORB782099-48, BYT-ORB782099-96 |
| Manufacturer | Biorbyt |
| Category | Kits/Assays |
| Application | ELISA |
| Species Reactivity | All |
| Product Description | This assay employs the competitive inhibition enzyme immunoassay technique. The microtiter plate provided in this kit has been pre-coated with Tyrosol(Tyr) protein. Standards or samples are added to the appropriate microtiter plate wells then with a ... |
| Concentration | 30000 ng/mL |
| Range | 468.75-30000 ng/mL |
| Sensitivity | 136.2 ng/mL |
| Samples | biological agents and other biological fluids |

Application Notes

Application Notes: standard: 30000 ng/mL. Test principle: This assay employs the competitive inhibition enzyme immunoassay technique. The microtiter plate provided in this kit has been pre-coated with Tyr protein. Standards or samples are added to the appropriate microtiter plate wells then with a biotin-conjugated antibody specific to Tyr. Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm 10nm. The concentration of Tyr in the samples is then determined by comparing the OD of the samples to the standard curve