

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### Anti-Metallothionein 3 Polyclonal Antibody TGM-TMAB-08763

|                            |                                                                                                                                                                                |
|----------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Article Name               | Anti-Metallothionein 3 Polyclonal Antibody                                                                                                                                     |
| Biozol Catalog Number      | TGM-TMAB-08763                                                                                                                                                                 |
| Supplier Catalog Number    | TMAB-08763                                                                                                                                                                     |
| Alternative Catalog Number | TGM-TMAB-08763-50UL,TGM-TMAB-08763-100UL,TGM-TMAB-08763-200UL                                                                                                                  |
| Manufacturer               | TargetMol                                                                                                                                                                      |
| Category                   | Biochemikalien                                                                                                                                                                 |
| Product Description        | Anti-Metallothionein 3 Polyclonal Antibody is a Rabbit antibody targeting Metallothionein 3. Anti-Metallothionein 3 Polyclonal Antibody can be used in IHC-P,IHC-Fr,IF,FCM.... |