

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### **Transferrin (TF, Beta-1 Metal-binding Globulin, Serotransferrin, Siderophilin) (Not for Export EU), Rabbit USB-T8199-01C**

|                            |                                                                                                                                                                                                                                                               |
|----------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Article Name               | Transferrin (TF, Beta-1 Metal-binding Globulin, Serotransferrin, Siderophilin) (Not for Export EU), Rabbit                                                                                                                                                    |
| Biozol Catalog Number      | USB-T8199-01C                                                                                                                                                                                                                                                 |
| Supplier Catalog Number    | T8199-01C                                                                                                                                                                                                                                                     |
| Alternative Catalog Number | USB-T8199-01C-1                                                                                                                                                                                                                                               |
| Manufacturer               | US Biological                                                                                                                                                                                                                                                 |
| Host                       | Rabbit                                                                                                                                                                                                                                                        |
| Category                   | Antikörper                                                                                                                                                                                                                                                    |
| Application                | ELISA                                                                                                                                                                                                                                                         |
| Immunogen                  | Purified transferrin isolated from pooled swine serum.                                                                                                                                                                                                        |
| Product Description        | Applications: Suitable for use in ELISA, Immunoelectrophoresis and Single and Double Radial Immunodiffusion. Other applications not tested. Recommended Dilution: Immunoelectrophoresis: 2ml plasma or equivalent against 120ml antiserum Double Radial Im... |
| Purity                     | Serum                                                                                                                                                                                                                                                         |
| Form                       | Supplied as a lyophilized powder in PBS, pH 7.2. Preservative free. Reconstitute with 1ml sterile dH2O.                                                                                                                                                       |