



Diagnostica Vertrieb GmbH, Oehleckerring 11-13

22419 Hamburg, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

QuantiQuik(TM) L-Lactic Acid Quick Test Strips BAS-QQLLAC10

Article Name	QuantiQuik(TM) L-Lactic Acid Quick Test Strips
Biozol Catalog Number	BAS-QQLLAC10
Supplier Catalog Number	QQLLAC10
Alternative Catalog Number	BAS-QQLLAC10
Manufacturer	BioAssay Systems
Category	Sonstiges
Product Description	For rapid, direct determination of L-lactic acid concentrations in food and beverage samples as well as biological samples. Semi-quantitative measurement between 0-360 mg/L (undiluted) L-lactic acid. Convenient. No expensive lab equipment needed. Sam
Storage	RT