



Diagnostica Vertrieb GmbH, Oehleckerring 11-13

22419 Hamburg, Germany

**Telephone:** +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## QuantiQuik(TM) Histamine Quick Test Strips BAS-QQHIST10

Article Name	QuantiQuik(TM) Histamine Quick Test Strips
Biozol Catalog Number	BAS-QQHIST10
Supplier Catalog Number	QQHIST10
Alternative Catalog Number	BAS-QQHIST10
Manufacturer	BioAssay Systems
Category	Sonstiges
Product Description	For rapid, direct determination of Histamine concentrations in fish, wine, and other samples. Key Features: Fast and sensitive. Use of 500 mg of fish or 250 $\mu$ L of sample. Semi-quantitative measurement between 0-200 ppm (undiluted) Histamine. Convenie
Storage	RT