

Diagnostica Vertrieb GmbH, Oehleckerring 11-13

22419 Hamburg, Germany

**Telephone:** +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## Anti-GAPDHS Purified, Clone: [Hs-8], Monoclonal EXB-11-301-C025

Article Name	Anti-GAPDHS Purified, Clone: [Hs-8], Monoclonal
Biozol Catalog Number	EXB-11-301-C025
Supplier Catalog Number	11-301-C025
Alternative Catalog Number	EXB-11-301-C025
Manufacturer	EXBIO
Category	Antikörper
Application	WB, FC, ICC
Species Reactivity	Porcine, Human
Immunogen	Freshly ejaculated human sperms were washed in PBS and extracted in 3% acetic acid, 10% glycerol, 30 mM benzaminidine. The acid extract was dialyzed against 0.2% acetic acid and subsequently used for immunization.
Product Description	GAPDHS (the sperm-specific glyceraldehyde phosphate dehydrogenase, also known as GAPD2, GAPDS, HSD-35, or GAPDH-2, is a glycolytic enzyme that plays an important role in carbohydrate metabolism. Like its somatic cell counterpart, this sperm-specific
Clonality	Monoclonal
Concentration	1 mg/ml
Clone Designation	[Hs-8]
Isotype	IgM
Buffer	Tris buffered saline (TBS), pH 8.0, 15 mM sodium azide

Storage	2°C to 8°C
Antibody Type	Monoclonal Antibody
Application Notes	Immunocytochemistry: Recommended dilution: 10 $\mu$ g/ml, membrane permeabilization (acetone) is essential. Flow cytometry: Recommended dilution: 3-12 $\mu$ g/ml. Intraacrosomal staining.