

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### Anti-ZAP70 FITC, Clone: [1E7.2], Monoclonal EXB-1F-665-T100

|                            |   |
|----------------------------|---|
| Article Name               | Anti-ZAP70 FITC, Clone: [1E7.2], Monoclonal   |
| Biozol Catalog Number      | EXB-1F-665-T100   |
| Supplier Catalog Number    | 1F-665-T100   |
| Alternative Catalog Number | EXB-1F-665-T100   |
| Manufacturer               | EXBIO   |
| Category                   | Antikörper  |
| Application                | FC  |
| Species Reactivity         | Human, Mouse  |
| Immunogen                  | A KLH-coupled peptide corresponding to amino acids 282-307 of human ZAP70   |
| Conjugation                | FITC  |
| Product Description        | The ZAP70 (zeta-associated protein of 70 kDa) tyrosine kinase was identified as a tyrosine phosphoprotein that associates with TCR zeta subunit and undergoes tyrosine phosphorylation following TCR stimulation. ZAP70 is a Syk family tyrosine kinase pr... |
| Clonality                  | Monoclonal  |
| Clone Designation          | [1E7.2]   |
| Isotype                    | Mouse IgG1 kappa  |
| Buffer                     | Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide   |

|                    |   |
|--------------------|---|
| Storage            | 2°C to 8°C  |
| Target             | ZAP70   |
| Antibody Type      | Monoclonal Antibody   |
| Application Dilute | Flow cytometry: The reagent is designed for analysis of human blood cells using 10 µl reagent / 100 µl of whole blood or 10 <sup>6</sup> cells in a suspension. The content of a vial (1 ml) is sufficient for 100 tests.                         |
| Application Notes  | Flow cytometry: The reagent is designed for analysis of human blood cells using 10 µl reagent / 100 µl of whole blood or 10 <sup>6</sup> cells in a suspension. The content of a vial (1 ml) is sufficient for 100 tests. Intracellular staining. |