

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Trident Membrane Protein Extraction Kit GTX-GTX16373

Article Name	Trident Membrane Protein Extraction Kit
Biozol Catalog Number	GTX-GTX16373
Supplier Catalog Number	GTX16373
Alternative Catalog Number	GTX16373-5,GTX16373-20
Manufacturer	GeneTex
Category	Kits/Assays
Application	WB, ELISA, IP
Application Dilute	1. Read the entire protocol carefully. 2. Before experiment, thaw buffer A or/and buffer B completely, invert the bottles a few times and keep them on ice. Place the filter cartridges in collection tubs, and pre-chill on ice. Pre-chill PBS on ice. <
Application Notes	1. Read the entire protocol carefully. 2. Before experiment, thaw buffer A or/and buffer B completely, invert the bottles a few times and keep them on ice. Place the filter cartridges in collection tubs, and pre-chill on ice. Pre-chill PBS on ice. <