



Diagnostica Vertrieb GmbH, Oehleckerring 11-13

22419 Hamburg, Germany

**Telephone:** +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## EnzyChrom(TM) NAD/NADH Assay Kit BAS-E2ND-100

Article Name	EnzyChrom(TM) NAD/NADH Assay Kit
Biozol Catalog Number	BAS-E2ND-100
Supplier Catalog Number	E2ND-100
Alternative Catalog Number	BAS-E2ND-100
Manufacturer	BioAssay Systems
Category	Kits/Assays
Product Description	For sensitive determination of NAD and NADH and evaluation of drug effects on NAD/NADH metabolism. Key Features: Sensitive and accurate. Detection limit of 0.05 $\mu\text{M}$ and linearity up to 10 $\mu\text{M}$ NAD+/NADH in 96-well plate assay. Convenient. The procedure