



Diagnostica Vertrieb GmbH, Oehleckerring 11-13

22419 Hamburg, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

EnzyChrom(TM) Acetaldehyde Assay Kit BAS-EACT-100

Article Name	EnzyChrom(TM) Acetaldehyde Assay Kit
Biozol Catalog Number	BAS-EACT-100
Supplier Catalog Number	EACT-100
Alternative Catalog Number	BAS-EACT-100
Manufacturer	BioAssay Systems
Category	Kits/Assays
Product Description	For quantitatve determination of acetaldehyde enzyme activity and evaluation of drug modulators. Key Features: Fast and sensitive. Linear detection range (20 μ L sample): 2 μ M to 2 mM acetaldehyde in 96-well plate assay. High-throughput.The procedure in