



Diagnostica Vertrieb GmbH, Oehleckerring 11-13

22419 Hamburg, Germany

**Telephone:** +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## EnzyChrom(TM) Adipolysis Assay Kit BAS-EAPL-200

Article Name	EnzyChrom(TM) Adipolysis Assay Kit
Biozol Catalog Number	BAS-EAPL-200
Supplier Catalog Number	EAPL-200
Alternative Catalog Number	BAS-EAPL-200
Manufacturer	BioAssay Systems
Category	Kits/Assays
Product Description	For quantitative assay of adipolysis and evaluation of drug effects on adipolysis. Key Features: Sensitive and accurate. Use as little as 10 $\mu L$ samples. Linear detection range in 96-well plate: 0.92 to 100 $\mu g/mL$ (10 to 1000 $\mu M$ ) glycerol for colorimet