



Diagnostica Vertrieb GmbH, Oehleckerring 11-13

22419 Hamburg, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

EnzyFluo(TM) ERK Phosphorylation Assay Kit BAS-EERK-100

Article Name	EnzyFluo(TM) ERK Phosphorylation Assay Kit
Biozol Catalog Number	BAS-EERK-100
Supplier Catalog Number	EERK-100
Alternative Catalog Number	BAS-EERK-100
Manufacturer	BioAssay Systems
Category	Kits/Assays
Product Description	For quantitative fluorescent immunoenzymatic assay of ERK1/2 phosphorylation status in cultured cells. Key Features: New and improved. Total assay time reduced from the standard 21 hours to 6.5 hours (hands-on time 2.5 hrs). Simple and convenient. Ce