



Diagnostica Vertrieb GmbH, Oehleckerring 11-13

22419 Hamburg, Germany

**Telephone:** +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## EnzyChrom(TM) Glycerol Assay Kit BAS-EGLY-200

Article Name	EnzyChrom(TM) Glycerol Assay Kit
Biozol Catalog Number	BAS-EGLY-200
Supplier Catalog Number	EGLY-200
Alternative Catalog Number	BAS-EGLY-200
Manufacturer	BioAssay Systems
Category	Kits/Assays
Product Description	For quantitative determination of glycerol and evaluation of drug effects on glycerol metabolism. Key Features: Sensitive and accurate. Use as little as 10 $\mu$ L samples. Linear detection range in 96-well plate: 10 to 1000 $\mu$ M (92 $\mu$ g/dL to 9.2 mg/dL) gly