



Diagnostica Vertrieb GmbH, Oehleckerring 11-13

22419 Hamburg, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

EnzyChrom(TM) Histamine Assay Kit BAS-EHIS-100

Article Name	EnzyChrom(TM) Histamine Assay Kit
Biozol Catalog Number	BAS-EHIS-100
Supplier Catalog Number	EHIS-100
Alternative Catalog Number	BAS-EHIS-100
Manufacturer	BioAssay Systems
Category	Kits/Assays
Product Description	For quantitative determination of histamine in food, beverage and agricultural products. Key Features: Fast and sensitive. Use of 20 μ L sample. Linear detection range 0.0025 - 1 mM (0.3 - 111 ppm) histamine in 96-well plate assay. Convenient. The pro