



Diagnostica Vertrieb GmbH, Oehleckerring 11-13

22419 Hamburg, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

EnzyChrom(TM) Pyruvate Assay Kit BAS-EPYR-100

Article Name	EnzyChrom(TM) Pyruvate Assay Kit
Biozol Catalog Number	BAS-EPYR-100
Supplier Catalog Number	EPYR-100
Alternative Catalog Number	BAS-EPYR-100
Manufacturer	BioAssay Systems
Category	Kits/Assays
Product Description	For quantitative determination of pyruvate and evaluation of drug effects on its metabolism. Key Features: Sensitive and accurate. Use as little as 10 μ L samples. Linear detection range in 96-well plate: 2 to 500 μ M (17 μ g/dL to 4.4 mg/dL) pyruvate f