



Diagnostica Vertrieb GmbH, Oehleckerring 11-13

22419 Hamburg, Germany

**Telephone:** +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## PhosphoBlocker(TM) Blocking Reagent CBL-AKR-104

Article Name	PhosphoBlocker(TM) Blocking Reagent
Biozol Catalog Number	CBL-AKR-104
Supplier Catalog Number	AKR-104
Alternative Catalog Number	CBL-AKR-104
Manufacturer	Cell Biolabs
Category	Sonstiges
Product Description	Most commercially available Western blot blockers, such as dry milk or serum, are sufficient to block unreactive sites on the membrane. However, they are not designed to preserve phosphoprotein antigens during blotting. Our PhosphoBLOCKER(TM) Blockin