



Diagnostica Vertrieb GmbH, Oehleckerring 11-13

22419 Hamburg, Germany

**Telephone:** +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## OxiSelect(TM) TBARS Assay Kit (MDA Quantitation) CBL-STA-330

Article Name	OxiSelect(TM) TBARS Assay Kit (MDA Quantitation)
Biozol Catalog Number	CBL-STA-330
Supplier Catalog Number	STA-330
Alternative Catalog Number	CBL-STA-330
Manufacturer	Cell Biolabs
Category	Kits/Assays
Application	Colorimetric/Fluorometric
Product Description	The TBARS (Thiobarbituric Acid Reactive Substances) assay is well-established for screening and monitoring lipid peroxidation. MDA forms a 1:2 adduct with thiobarbituric acid, the MDA-TBA adduct can then be measured. Our OxiSelect(TM) TBARS Assay Kit