

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

OxiSelect(TM) TBARS Assay Kit (MDA Quantitation) CBL-STA-330-5

| | |
|----------------------------|---|
| Article Name | OxiSelect(TM) TBARS Assay Kit (MDA Quantitation) |
| Biozol Catalog Number | CBL-STA-330-5 |
| Supplier Catalog Number | STA-330-5 |
| Alternative Catalog Number | CBL-STA-330-5 |
| Manufacturer | Cell Biolabs |
| Category | Kits/Assays |
| Application | Colorimetric/Fluorometric |
| Product Description | The TBARS (Thiobarbituric Acid Reactive Substances) assay is well-established for screening and monitoring lipid peroxidation. MDA forms a 1:2 adduct with thiobarbituric acid, the MDA-TBA adduct can then be measured. Our OxiSelect(TM) TBARS Assay Kit... |