



Diagnostica Vertrieb GmbH, Oehleckerring 11-13

22419 Hamburg, Germany

**Telephone:** +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## Pan-Ras Activation Assay CBL-STA-400

Article Name	Pan-Ras Activation Assay
Biozol Catalog Number	CBL-STA-400
Supplier Catalog Number	STA-400
Alternative Catalog Number	CBL-STA-400
Manufacturer	Cell Biolabs
Category	Kits/Assays
Application	WB
Product Description	Our Ras Activation Assays use visible agarose beads to selectively precipitate the active form of specific Ras protein of interest. The precipitated small GTPase is then detected by Western blot using a target-specific antibody included in the kit. A