



Diagnostica Vertrieb GmbH, Oehleckerring 11-13

22419 Hamburg, Germany

**Telephone:** +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## Arf6 Activation Assay CBL-STA-407-6

Article Name	Arf6 Activation Assay
Biozol Catalog Number	CBL-STA-407-6
Supplier Catalog Number	STA-407-6
Alternative Catalog Number	CBL-STA-407-6
Manufacturer	Cell Biolabs
Category	Kits/Assays
Application	WB
Product Description	Our Arf Activation Assays use visible agarose beads to selectively precipitate the active form of Arf1 or Arf 6. The precipitated small GTPase is then detected by Western blot using an Arf1- or Arf6-specific antibody included in the kit