

Diagnostica Vertrieb GmbH, Oehleckerring 11-13

22419 Hamburg, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Anti-PsaO | PSI-O subunit of photosystem I, Rabbit, Polyclonal AGR-AS04-050

Article Name	Anti-PsaO PSI-O subunit of photosystem I, Rabbit, Polyclonal
Biozol Catalog Number	AGR-AS04-050
Supplier Catalog Number	AS04-050
Alternative Catalog Number	AGR-AS04-050
Manufacturer	Agrisera
Host	Rabbit
Category	Antikörper
Application	WB
Species Reactivity	A. thaliana, Plant
Immunogen	Fusion protein between DHFR and the mature part of PSI-O of Arabidopsis thaliana UniProt: Q949Q5TAIR: At1g08380
Product Description	PsaO - photosystem I subunit O. The mature PsaO is a 10-kDa protein with two transmembrane helices. It has no counterpart in Photosystem I of cyanobacteria but seems to be present in higher plants, in mosses, and in green algae
Clonality	Polyclonal
Molecular Weight	15.1 10.1 kDa
NCBI	837358
UniProt	Q8L7C3
Purity	Total IgG. Protein G purified in PBS pH 7.4.

Form	Lyophilized
Antibody Type	Polyclonal Antibody
Application Dilute	1 : 2000 - 5000 (WB)
Application Notes	For reconstitution add 100 µl of sterile water