



Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

InjectMan 4, micromanipulator, power plug EU EPE-5192000019

Article Name	InjectMan 4, micromanipulator, power plug EU
Biozol Catalog Number	EPE-5192000019
Supplier Catalog Number	5192000019
Alternative Catalog Number	EPE-5192000019
Manufacturer	BIOZOL
Category	Zellmanipulation
Product Description	InjectMan 4 Minimal lesion hole size for high cell survival rates are key in your DNA/RNA, protein or fluorescent dye microinjection experiment? Or you want to speed up your serial adherent cell injections to ~20-30 cells/min? The InjectMan 4 elec